



Vertrieb durch:

Immundiagnostik AG

Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

www.immundiagnostik.com

E-Mail: info@immundiagnostik.com



IM-CCA030-140718E

FOR RESEARCH USE ONLY

Catalog #CCA-030

MC Reagent Biochemical Assay Series

Calcium Assay Kit

O-Cresolphthalein (OCPC) Chromogenic method

Biochemical Significance and Test Summary

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones. The calcium level in the extracellular space is in dynamic equilibrium with bone calcium. Calcium ions affect the contractility of heart and skeletal muscles and are essential for the function of the nervous system. Calcium ions play an important role in blood clotting and bone minimization. In plasma, calcium is bound to a considerable extent to proteins, 10 % is in the form of inorganic complexes and 50 % is present as free ion species. The calcium homeostasis regulated by the parathyroid hormone (PTH), calcitriol (CT), and calcitonin. A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption. Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism.

This product is a direct colorimetric assay kit without deproteinization of the sample. Calcium with o-cresolphthale (as chelator) at strong alkaline pH, yields a violet colored complex. The intensity of the color formed is proportional to the calcium concentration in the sample.

1. Kit contents (for 200 tests)

R-1	Buffer	1 x 38 mL	Ready to use
R-2	Chelate color	1 x 10 mL	Ready to use
STD	10 mg/dL Ca standard	1 x 1.0 mL	Ready to use

*Storage conditions: Store at 2-8°C. **Don't freeze.**

*Expiration: 1 year. After the vials are opened, the kit should be used in one month.

*Measuring range: 0.2-30 mg/dL

2. Materials required but not provided

- (1) Distilled water
- (2) Micropipettors and pipette tips
- (3) Clear flat-bottom 96-well plate
- (4) Microplate reader with 570 nm capability

3. Assay preparation

Bring all reagents to room temperature before use.

4. Sample preparation

Serum/Plasma: Insoluble substances in serum and plasma samples should be removed by filtration or centrifugation. EDTA-plasma cannot be used.

Urine (24 hour pooled urine): Dilute the sample 1/2 in distilled water. Add 6M HCl to the diluted sample and adjust pH 2.0-3.0 (e.g. 5-10µL of 6M HCl/ 1mL of lysate). Centrifuge at 6,000 rpm for 15 minutes. Collect the supernatant and use it for assay.

Biological fluid: Add 6M HCl to the sample and adjust pH 2.0-3.0 (e.g. 5-10µL of 6M HCl/ 1mL of lysate). Centrifuge at 6,000 rpm for 15 minutes. Collect the supernatant and use it for assay.

Tissue: Add 3% TCA solution, vortex for 1 minute and incubate for 30 minutes at 4-8°C. Centrifuge at 6,000 rpm for 15 minutes. Collect the supernatant and use it for assay.

Note: Sample pH should be between pH 2 and pH 8.

5. Assay protocol

- (1) Add 5 µL of Distilled water(Blank)/STD(Standard)/Sample into each well.

- (2) Add 190 μ L of R-1 to each well and incubate for 5 minutes at room temperature.
- (3) Add 50 μ L of R-2 to each well and incubate for 5 minutes at room temperature.
- (4) Read the absorbance at 570 nm (560-590 nm). ----- OD

6. Calculations

$$\Delta OD_{\text{Standard}} = OD_{\text{Standard}} - OD_{\text{Blank}}, \Delta OD_{\text{Sample}} = OD_{\text{Sample}} - OD_{\text{Blank}}$$

$$\text{Calcium (mg/dL)} = \Delta OD_{\text{Sample}} / \Delta OD_{\text{Standard}} \times 10$$

$$\text{Calcium (mM)} = \Delta OD_{\text{Sample}} / \Delta OD_{\text{Standard}} \times 2.495$$

(Assay example)

	OD (570 nm)	Δ OD	Calcium (mg/L)
DW (Blank)	0.677	-	-
Standard	1.116	0.439	-
Sample	0.976	0.299	6.81

$$\Delta OD_{\text{Standard}} = 1.116 - 0.677 = 0.439$$

$$\Delta OD_{\text{Sample}} = 0.976 - 0.677 = 0.299$$

$$\text{Calcium}_{\text{Sample}} = \Delta OD_{\text{Sample}} / \Delta OD_{\text{Standard}} \times 10 = (0.299 / 0.439) \times 10 = 6.81 \text{ (mg/dL)}$$

$$\text{Calcium}_{\text{Sample}} (\mu\text{M}) = \Delta OD_{\text{Sample}} / \Delta OD_{\text{Standard}} \times 2.495 = (0.299 / 0.439) \times 2.495 = 1.70 \text{ (mM)}$$

7. Interferences

EDTA inhibits calcium to chromogenic system. The test is not affected by presence of bilirubin-F and bilirubin-C up to 40 mg/dL and chyle up to 1,000 FTU.

8. Quality Control

Use of control sera is recommended to monitor the quality of assay results.

9. Reference

- (1) J. W. Ferguson, J. J. Richard, J. W. O'laughlin and C. V. Banks: Simultaneous Spectrophotometric Determination of Calcium and Magnesium with Chlorophosphonazo-III, *Anal. Chem*, 36, 796.2 (1962).
- (2) Harold. V, Connerty. MD, Anglis. R, Briggs. BS: Determination of serum calcium by means of orthocresolphthalein complexone. *Am J Clin Pathol*, 45, p290-296 (1966).
- (3) Fujita. T, Noguchi. K, Terashima. I : Apoplastic mesophyll signals induce rapid stomatal responses to CO₂ in *Commelina communis*, *New Phytol*, 199(2), p395-406 (2012).

10. Technical support & troubleshooting

- (1) Unstablensness of incubation temperature may result in unstable results.
- (2) Use disposable test tube and glassware washed with 1M HNO₃ or 1M HCl, and rinse with distilled water.
- (3) Accuracy to the microliter is important to obtain good results. Ensure maximum precision when pipetting.
- (4) Temperature for the chromogenic reaction may affect the optical density. It may be necessary to adjust the reaction time depending on the room temperature.
- (5) High concentration of proteins or lipid in cell lysate or in tissue extract may affect the observed value. Please remove them by ultrafiltration or centrifugation.
- (6) 24 hour pooled urine: Dilute the sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor) for assay sample.



**Japan Institute for the Control of Aging (JaICA),
Nikken SEIL Co., Ltd.**

Haruoka 710-1, Fukuroi, Shizuoka 437-0122, Japan

TEL: +81-538-49-0125 FAX: +81-538-49-1267

www.jaica.com/e/ E-mail: biotech@jaica.com

(Supported languages: Japanese and English)